

KOALA DAQ Operation Manual

1. Login method & passwords

- a. **Control room monitor server:**
 - i. **User name:** koala
 - ii. **Password:** only4koala
- b. **DAQ server:** `ssh -Y koala-daq`
- c. **Online server:** `ssh -Y koala-online`

2. Environment Setup:

First log into one of the control room PC, then:

- a. Setup the running environment for DAQ:
 - i. In **DAQ server terminal**, type command: **ems_setup**
- b. *Open DAQ GUI:*
 - i. In **DAQ server terminal**, type command: **ems_control**
- c. *Open Online program:*
 - i. In **Online server terminal**, type command: **ems_online**

3. Data taking Operation (all in *ems_control* GUI):

In each data taking, click the following buttons sequentially:

- a. **INIT:**
 - i. to initialize the system
 - ii. “**Init complete**” on the GUI log-window indicates successful initialization
 - iii. **NOTE: START** only after a successful **INIT**
- b. **START:**
 - i. to start the data taking run
 - ii. A new data file will be created, the file name is shown under “**file name**”
 - iii. **NOTE: START** must follow a successful **INIT**
- c. **STOP:**
 - i. to stop the current data taking run
 - ii. “**Stop of \$run_nr complete**” on the GUI log-window indicates successful stop
 - iii. The stopping process may take a while to release the resources (~60s)
 - iv. **NOTE: INIT** the system only after a successful **STOP** (i.e. wait ~60 seconds before re-INIT)

4. In case of beam momentum changing, also change the associated filename prefix

- a. In DAQ server terminal, type: **ems_ch_mom \$new_momentum_value**
- b. In DAQ GUI, click menu: “**File**” --> “**Restor Setup**”
- c. Then, the normal data taking operation: **INIT --> START --> STOP**

KOALA Logging Check-List

1. Date & Time

2. Beam condition

- Momentum (ask COSY crew)
- Intensity (protons in COSY ring)
- Scaler values (from WebCam, ssh to koala-daq, run kamerka on terminal)

3. Target condition

- Remote desktop: (the curves should be stable)
- Gas System status: ***koalatarget.ikp.kfa-juelich.de***

AG Khroukaz		MCT-K target	Impressum
Gas System			
Nozzle Baratron	15.77	bar	
Gas Flow	2.02	l/min	
Nozzle Temperature	35.00	K	
Nozzle Diameter	35.77	µm	
Cluster Source			
Skimmer Chamber	9.22e-2	mbar	
Collimator Chamber	6.10e-4	mbar	
Interaction Chamber			
Interaction Chamber	0.00e+0	mbar	
Beam Dump			
Beam Dump 1	7.87e-6	mbar	
Beam Dump 2	1.24e-6	mbar	
Beam Dump 3	2.34e-5	mbar	

4. Detector condition

- Sensor temperature (command line: ***lakeshore336***)
- High voltage
 - Recoil (command line: ***mpod***)
 - Fwd (command line: ***telnet 134.94.192.85 1527***)
- Vacuum level (IPcamera: ***<http://134.94.192.148>, name: admin, passwd: only4koala***)
- Oscilloscope remote access: ***xfreerdp -sec-nla /u:Tek_Local_Admin /p:only4koala /v:134.94.192.152:3389***

5. Data taking

- Start & Stop time
- Run number
- File name
- Checking all the conditions above in the list